

### Hand-held LLLT

THOR Laser hand-held devices deliver focused beams of Low-Level Light Therapy (LLLT), also known as Photobiomodulation (PBM) Therapy, to reduce pain and inflammation, repair tissue, and promote healing at the cellular level.

PBM Therapy treatments take only a few minutes and can be applied two or more times a week. Applications include:

- Sprains and strains
- Osteoarthritis
- Post-operative pain
- Pain and inflammation
- Injury prevention
- Injury and post-surgical rehab
- Localized swelling
- Enhancing performance

### How It Works

When cells undergo stress from causes such as inadequate diet, smoking, or injury, the result is an increased amount of free radicals that inhibit proper cell function including the production of energy. PBM Therapy works by displacing free radicals to reduce oxidative stress, increase circulation, and increase cellular energy production. It's a safe, effective, and drug-free modality with no known side effects.



# THOR Laser | Specifications

## TECHNICAL DATA

LX2 table top control unit

Large healing / anti-inflammatory 69 LED cluster

Analgesic 810nm 1Watt laser cluster (5 x 200mW)

Trigger Point 810nm 200mW laser point locator

2 x Laser safety glasses

## WARRANTY

All THOR Laser devices with a 5-year limited warranty on all parts and labor, and access to the THOR treatment library.

## PRICING

Please call Sales at 214.363.3111 for all pricing information.

