The CryoNiQ Cryochamber is a single-person chamber used for individual sessions of Whole Body Cryotherapy. During a session, clients wear minimal protective clothing (gloves, socks, and briefs). The temperature is regulated by the automated system but can also be adjusted within a range of -166°F to -320°F by the operating technician via keypad or optional touchscreen control panel.

The unit is equipped with an electric lift used to adjust a client's height in the chamber so that the body is submerged in the hypercooled air up to the neck while the head remains above. This enables the client to be in continuous communication with the technician, which improves the client's overall comfort.

For added safety, the CryoNiQ Cryochamber comes standard with a pair of non-locking pivoting doors; one or both can be opened at any time. The unit also features a heated ventilating system to minimize moisture residues and to vent out used nitrogen gasiform air. The geometric design of the inner cabin optimizes convection of cryotherapy temperatures for improved session results.

- Temperature Range: -166°F to -320°F
- Continuous Flow Technology for stable temperature throughout the session
- Electric lift with push button control
- Unique dual-door design for easier entering and exiting
- Customizable with multiple exterior and interior color options (Icy White shown)
- Session time: 2-3 min
- Nitrogen usage during session: avg of 6 liters per session
- E-Tablet Control Panel and viewing screen
- Automatic standby 4 cooling





## Product Specifications CryoNiQ Cryochamber

#### **TECHNICAL DATA**

- Cryogenic gas: Liquid Nitrogen
- A/C Power: 220V, single-phase electricity
- Width (doors open): 84 in.
- Width (doors closed): 37 in.
- Length: 59 in.
- · Height: 95 in.
- Diameter of chamber: 37 in.
- Dimension of the hydrogen tank: 28"L x 29"W x 55"H

#### **WARRANTY**

All CryoNiQ Cryochambers come with a 2-year limited warranty on all parts and labor.

### **PRICING**

Please call Sales at 214.363.3111 for all pricing information.

# 180 **Client Viewing Screen Operator Viewing Screen FINISH OPTIONS**

Exterior	Interior
Icy White Light Grey Cashmere Dark Grey Black	Red Black Grey Navy Blue Blue Refl.
Dark Teal Carmine Red  Standard Exterior: Light Grey.	Orange Refl.  Standard Interior:  Black.  All other options, \$1000

All other options, \$1000.



Cold therapy products and equipment have not been tested or approved by the FDA or any other government agencies for the treatment of any illness or disease. Use at your own risk.

